

PDGFR α Mouse Monoclonal Antibody(4G11)

Catalog No: #44208

Package Size: #44208-1 50ul #44208-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	PDGFR α Mouse Monoclonal Antibody(4G11)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	4G11
Isotype	IgG1
Applications	IHC
Species Reactivity	Hu Rt Ms
Specificity	PDGFR α Mouse Monoclonal Antibody(4G11) detects endogenous levels of PDGFRA
Immunogen Type	Peptide
Immunogen Description	Synthetic Peptide
Conjugates	Unconjugated
Target Name	PDGFR α
Other Names	PDGFRA
Accession No.	Swiss-Prot:P16234Gene ID:5156
Calculated MW	180kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Immunohistochemistry: 1:100-1:200

Images



Immunohistochemical analysis of paraffin-embedded Rat Skin Tissue using PDGFR α Mouse mAb diluted at 1:200.



Immunohistochemical analysis of paraffin-embedded Human Skin Tissue using PDGFR a Mouse mAb diluted at 1:200.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.