

JAK2 Mouse Monoclonal Antibody(9B4)

Catalog No: #44248

Package Size: #44248-1 50ul #44248-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	JAK2 Mouse Monoclonal Antibody(9B4)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	9B4
Isotype	IgG1
Applications	IHC
Species Reactivity	Hu Ms Rt
Specificity	The antibody detects endogenous JAK2 protein
Immunogen Type	Peptide
Immunogen Description	Synthetic Peptide of JAK2
Conjugates	Unconjugated
Target Name	JAK2
Other Names	Tyrosine-protein kinase JAK2 (EC 2.7.10.2) (Janus kinase 2) (JAK-2)
Accession No.	Swiss-Prot:O60674Gene ID:3717
Calculated MW	120kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Immunohistochemistry: 1:50-1:200

Images



Immunohistochemical analysis of paraffin-embedded Human Breast Carcinoma Tissue using JAK2 Mouse mAb diluted at 1:200



Immunohistochemical analysis of paraffin-embedded Human Lung Carcinoma Tissue using JAK2 Mouse mAb diluted at 1:200

Note: This product is for in vitro research use only and is not intended for use in humans or animals.