

P38 Mouse Monoclonal Antibody(1G1)

Catalog No: #44249

Package Size: #44249-1 50ul #44249-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	P38 Mouse Monoclonal Antibody(1G1)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	1G1
Isotype	IgG1
Applications	IHC
Species Reactivity	Hu Ms Rt
Specificity	The antibody detects endogenous P38 protein
Immunogen Type	Peptide
Immunogen Description	Synthetic Peptide of P38
Conjugates	Unconjugated
Target Name	P38
Other Names	Mitogen-activated protein kinase 14 (MAP kinase 14) (MAPK 14) (EC 2.7.11.24) (Cytokine suppressive anti-inflammatory drug-binding protein) (CSAID-binding protein) (CSBP) (MAP kinase MXI2) (MAX-interacting protein 2) (Mitogen-activated protein kinase p38 alpha) (MAP kinase p38 alpha) (Stress-activate
Accession No.	Swiss-Prot:Q16539Gene ID:1432
Calculated MW	42kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Immunohistochemistry: 1:50-1:200

Images



Immunohistochemical analysis of paraffin-embedded Human Lung Carcinoma Tissue using P38 Mouse mAb diluted at 1:200

Note: This product is for in vitro research use only and is not intended for use in humans or animals.