

NBR1 Mouse Monoclonal Antibody(7C3)

Catalog No: #44253

Package Size: #44253-1 50ul #44253-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	NBR1 Mouse Monoclonal Antibody(7C3)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	7C3
Isotype	IgG1
Applications	IHC
Species Reactivity	Hu Ms Rt
Specificity	The antibody detects endogenous NBR1 protein
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant Protein of NBR1
Conjugates	Unconjugated
Target Name	NBR1
Other Names	Next to BRCA1 gene 1 protein (Cell migration-inducing gene 19 protein) (Membrane component chromosome 17 surface marker 2) (Neighbor of BRCA1 gene 1 protein) (Protein 1A1-3B)
Accession No.	Swiss-Prot:Q14596Gene ID:4077
Calculated MW	120kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Immunohistochemistry: 1:50-1:200

Images



Immunohistochemical analysis of paraffin-embedded Human Skeletal Muscle Tissue using NBR1 Mouse mAb diluted at 1:200.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.