

ATG5 Mouse Monoclonal Antibody(8B4)

Catalog No: #44255

Package Size: #44255-1 50ul #44255-2 100ul


Description

Product Name	ATG5 Mouse Monoclonal Antibody(8B4)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	8B4
Isotype	IgG1
Applications	IHC
Species Reactivity	Hu Ms Rt
Specificity	The antibody detects endogenous ATG5 protein
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant Protein of ATG5
Conjugates	Unconjugated
Target Name	ATG5
Other Names	Autophagy protein 5 (APG5-like) (Apoptosis-specific protein)
Accession No.	Swiss-Prot:Q9H1Y0Gene ID:9474
Calculated MW	55kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Immunohistochemistry: 1:50-1:200

Images



Immunohistochemical analysis of paraffin-embedded Human Ovarian Carcinoma Tissue using ATG5 Mouse mAb diluted at 1:200.



Immunohistochemical analysis of paraffin-embedded Human Breast Carcinoma Tissue using ATG5 Mouse mAb diluted at 1:200.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.