

OSGI2 rabbit pAb

Catalog No: #44256

Package Size: #44256-1 50ul #44256-2 100ul


Description

Product Name	OSGI2 rabbit pAb
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Applications	WB
Species Reactivity	Human; Mouse;Rat
Specificity	This antibody detects endogenous levels of OSGI2 at Human
Immunogen Description	Synthesized peptide derived from human OSGI2
Conjugates	Unconjugated
Accession No.	Swiss Prot:Q9Y236GeneID:734
Calculated MW	56kD
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

WB 1oO 500-2000

Images



Note: This product is for in vitro research use only and is not intended for use in humans or animals.