

PR Polyclonal Antibody

Catalog No: #44262

Package Size: #44262-1 50ul #44262-2 100ul

Description

Product Name	PR Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB;IHC;IF/ICC;ELISA
Species Reactivity	Human; Mouse;Rat
Specificity	PR Polyclonal Antibody detects endogenous levels of PR protein.
Immunogen Description	The antiserum was produced against synthesized peptide derived from human Progesterone Receptor.
Conjugates	Unconjugated
Accession No.	Swiss Prot:P06401GenelD:5241
Calculated MW	99kD
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.


Immunofluorescence: 1/200 - 1/1000.

ELISA: 1/5000.


Images



Western Blot analysis of various cells using PR Polyclonal Antibody



Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue, using Progesterone Receptor Antibody. The picture on the right is blocked with the synthesized peptide.



Immunofluorescence analysis of A549 cells, using Progesterone Receptor Antibody.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.