

ATG7 Mouse Monoclonal Antibody(3D6)

Catalog No: #44265

Package Size: #44265-1 50ul #44265-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ATG7 Mouse Monoclonal Antibody(3D6)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	3D6
Isotype	IgG1
Applications	IHC
Species Reactivity	Hu Ms Rt
Specificity	The antibody detects endogenous ATG7 protein
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant Protein of ATG7
Conjugates	Unconjugated
Target Name	ATG7
Other Names	Ubiquitin-like modifier-activating enzyme ATG7 (ATG12-activating enzyme E1 ATG7) (Autophagy-related protein 7) (APG7-like) (hAGP7) (Ubiquitin-activating enzyme E1-like protein)
Accession No.	Swiss-Prot:O95352Gene ID:10533
Calculated MW	78kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Immunohistochemistry: 1:50-1:200

Images



Immunohistochemical analysis of paraffin-embedded Human Breast Carcinoma Tissue using ATG7 Mouse mAb diluted at 1:200.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.