

PHF5A rabbit pAb

Catalog No: #44277

Package Size: #44277-1 50ul #44277-2 100ul

Description

| | |
|-----------------------|---|
| Product Name | PHF5A rabbit pAb |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Isotype | IgG |
| Purification | The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen. |
| Applications | WB;ELISA;IHC |
| Species Reactivity | Human; Mouse;Rat |
| Specificity | This antibody detects endogenous levels of PHF5A at Human/Mouse/Rat |
| Immunogen Description | Synthesized peptide derived from human PHF5A |
| Conjugates | Unconjugated |
| Accession No. | Swiss-Prot#: Q7RTV0 NCBI Gene#: 84844 |
| Calculated MW | 12kD |
| Concentration | 1 mg/ml |
| Formulation | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Storage | Store at -20°C/1 year |

Application Details

WB 1:500-2000

IHC-p 1:50-300

ELISA 2000-20000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.