

RASGEF1B Antibody

Catalog No: #46172

Package Size: #46172-1 50ul #46172-2 100ul



SAB[®]
Signalway Antibody

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	RASGEF1B Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	RASGEF1B Antibody detects endogenous levels of total RASGEF1B
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human RASGEF1B
Conjugates	Unconjugated
Target Name	RASGEF1B
Other Names	GPIG4B£B~Ras-GEF domain-containing family member 1B
Accession No.	Swiss-Prot#: Q0VAM2NCBI Gene ID:153020
Calculated MW	55 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Western Blot analysis of Hela lysates, using primary antibody at 1:1000 dilution.

Published Papers

Menglong Wang;Tongrui Wang;Yang Liu;Lurong Zhou;Yuanping Yin;Feng Gu et al., Identification and study of mood-related biomarkers and potential molecular mechanisms in type 2 diabetes mellitus., , (2025)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.