

TGFβ RIII Polyclonal Antibody

Catalog No: #46732

Package Size: #46732-1 50ul #46732-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	TGFβ RIII Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC
Species Reactivity	Hu
Specificity	TGFβ RIII Polyclonal Antibody detects endogenous levels of TGFβ RIII
Immunogen Type	protein
Immunogen Description	Synthesized peptide derived from TGFβ RIII at AA range: 731-780
Conjugates	Unconjugated
Target Name	TGFβ RIII
Other Names	transforming growth factor, beta receptor III
Accession No.	Swiss-Prot:Q03167NCBI Gene ID:7049
Calculated MW	93 kDa
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C


Application Details

Western blotting: 1B£B1500-1B£B12000

Images



Western Blot analysis of mouse-kidney mouse-heart mouse-brain mouse-lung cells using TGFβ RIII Polyclonal Antibody diluted at 1:1000. Secondary antibody was diluted at 1:20000



Immunohistochemical analysis of paraffin-embedded human-placenta, antibody was diluted at 1:200

Note: This product is for in vitro research use only and is not intended for use in humans or animals.