

CYFIP2 Polyclonal Antibody

Catalog No: #46735

Package Size: #46735-1 50ul #46735-2 100ul

Description


Product Name	CYFIP2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC
Species Reactivity	Hu Ms Rt
Specificity	CYFIP2 Polyclonal Antibody detects endogenous levels of CYFIP2
Immunogen Type	protein
Immunogen Description	Synthesized peptide derived from CYFIP2 at AA range: 1171-1220
Conjugates	Unconjugated
Target Name	CYFIP2
Other Names	Cytoplasmic FMR1-interacting protein 2 (p53-inducible protein 121)
Accession No.	Swiss-Prot:Q96F07NCBI Gene ID:26999
Calculated MW	195 kDa
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C

Application Details


Western blotting: 1B£B1500-1B£B12000

Immunohistochemistry: IHC: 1:100-1:300


Images




Western Blot analysis of HEPG2 cells using CYFIP2 Polyclonal Antibody diluted at 1:500. Secondary antibody was diluted at 1:20000



Immunohistochemical analysis of paraffin-embedded human-colon, antibody was diluted at 1:200



Immunohistochemical analysis of paraffin-embedded human-colon, antibody was diluted at 1:200



Immunohistochemical analysis of paraffin-embedded human-tonsil, antibody was diluted at 1:200

Note: This product is for in vitro research use only and is not intended for use in humans or animals.