

# ERLIN1/2 Polyclonal Antibody

Catalog No: #46738

Package Size: #46738-1 50ul #46738-2 100ul

## Description


|                       |                                                                                                                       |
|-----------------------|-----------------------------------------------------------------------------------------------------------------------|
| Product Name          | ERLIN1/2 Polyclonal Antibody                                                                                          |
| Host Species          | Rabbit                                                                                                                |
| Clonality             | Polyclonal                                                                                                            |
| Purification          | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Applications          | WB IHC                                                                                                                |
| Species Reactivity    | Hu Ms Rt                                                                                                              |
| Specificity           | ERLIN1/2 Polyclonal Antibody detects endogenous levels of ERLIN1/2                                                    |
| Immunogen Description | Synthesized peptide derived from ERLIN1/2 at AA range: 31-80                                                          |
| Conjugates            | Unconjugated                                                                                                          |
| Target Name           | ERLIN1/2                                                                                                              |
| Accession No.         | Swiss-Prot:O75477/O94905NCBI Gene ID:10613/11160                                                                      |
| Calculated MW         | 82 kDa                                                                                                                |
| Concentration         | 1 mg/ml                                                                                                               |
| Formulation           | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                               |
| Storage               | Store at -20°C                                                                                                        |

## Application Details


Western blotting: 1µg/500-1µg/2000

Immunohistochemistry: IHC: 1:100-1:300


## Images



Western Blot analysis of HEPG2 3T3 HeLa cells using ERLIN1/2 Polyclonal Antibody diluted at 1:800. Secondary antibody was diluted at 1:20000



Immunohistochemical analysis of paraffin-embedded human-liver, antibody was diluted at 1:200



Immunohistochemical analysis of paraffin-embedded human-liver, antibody was diluted at 1:200

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.