

Perforin 1 Polyclonal Antibody

Catalog No: #46739

Package Size: #46739-1 50ul #46739-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	Perforin 1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC
Species Reactivity	Hu
Specificity	Perforin 1 Polyclonal Antibody detects endogenous levels of Perforin 1
Immunogen Type	protein
Immunogen Description	Synthesized peptide derived from Perforin 1 at AA range: 451-500
Conjugates	Unconjugated
Target Name	Perforin 1
Other Names	Perforin-1 (P1) (Cytolysin) (Lymphocyte pore-forming protein) (PFP)
Accession No.	Swiss-Prot:P14222NCBI Gene ID:5551
Calculated MW	61 kDa
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C

Application Details


Western blotting: 1B£B1500-1B£B12000

Immunohistochemistry: IHC: 1:100-1:300


Images




Western Blot analysis of HeLa cells using Perforin 1 Polyclonal Antibody diluted at 1:500. Secondary antibody was diluted at 1:20000



Immunohistochemical analysis of paraffin-embedded human-kidney, antibody was diluted at 1:200



Immunohistochemical analysis of paraffin-embedded human-kidney, antibody was diluted at 1:200



Immunohistochemical analysis of paraffin-embedded human-spleen, antibody was diluted at 1:200

Note: This product is for in vitro research use only and is not intended for use in humans or animals.