

Imp3 Polyclonal Antibody

Catalog No: #46748

Package Size: #46748-1 50ul #46748-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	Imp3 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Applications	WB IHC
Species Reactivity	Hu Ms Rt
Specificity	The antibody detects endogenous Imp3 protein
Immunogen Type	protein
Immunogen Description	Synthetic peptide from human protein at AA range: 135-184
Conjugates	Unconjugated
Target Name	Imp3
Other Names	IMP3 C15orf12 MRPS4
Accession No.	Swiss-Prot:Q9NV31NCBI Gene ID:55272
Calculated MW	25 kDa
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C

Application Details


Western blotting: 1B£B1500-1B£B110000

Immunohistochemistry: IHC: 1:100-1:300


Images



Immunohistochemical analysis of paraffin-embedded human-kidney, antibody was diluted at 1:200



Immunohistochemical analysis of paraffin-embedded human-kidney, antibody was diluted at 1:200



Immunohistochemical analysis of paraffin-embedded human-lung-cancer, antibody was diluted at 1:200

Published Papers

el et., A novel tumour enhancer function of Insulin-like growth factor II mRNA-binding protein 3 in colorectal cancerInCell Death DisOn2023 Apr 6byDavide Di Fusco?1,?Antonio Di Grazia et al..PMID:?37024466, , (2023)

[PMID:37024466](#)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.