

## COL17A1 Polyclonal Antibody

Catalog No: #46757

Package Size: #46757-1 50ul #46757-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description


|                       |                                                                                                           |
|-----------------------|-----------------------------------------------------------------------------------------------------------|
| Product Name          | COL17A1 Polyclonal Antibody                                                                               |
| Host Species          | Rabbit                                                                                                    |
| Clonality             | Polyclonal                                                                                                |
| Purification          | The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen. |
| Applications          | WB IHC                                                                                                    |
| Species Reactivity    | Hu                                                                                                        |
| Specificity           | The antibody detects endogenous COL17A1 protein                                                           |
| Immunogen Type        | protein                                                                                                   |
| Immunogen Description | Synthetic peptide from human protein at AA range: 481-530                                                 |
| Conjugates            | Unconjugated                                                                                              |
| Target Name           | COL17A1                                                                                                   |
| Other Names           | COL17A1 BP180 BPAG2                                                                                       |
| Accession No.         | Swiss-Prot:Q9UMD9NCBI Gene ID:1308                                                                        |
| Calculated MW         | 150 kDa                                                                                                   |
| Concentration         | 1 mg/ml                                                                                                   |
| Formulation           | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                   |
| Storage               | Store at -20°C                                                                                            |

## Application Details


Western blotting: 1B£B1500-1B£B110000

Immunohistochemistry: IHC: 1:100-1:300

## Images



Immunohistochemical analysis of paraffin-embedded human-colon, antibody was diluted at 1:200



Immunohistochemical analysis of paraffin-embedded human-skin, antibody was diluted at 1:200

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.