

Tcl1 Polyclonal Antibody

Catalog No: #46760

Package Size: #46760-1 50ul #46760-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Tcl1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Applications	WB
Species Reactivity	Hu
Specificity	The antibody detects endogenous Tcl1
Immunogen Type	protein
Immunogen Description	Synthetic peptide from human protein at AA range: 30-110
Conjugates	Unconjugated
Target Name	Tcl1
Other Names	T-cell leukemia/lymphoma protein 1A (Oncogene TCL-1) (Oncogene TCL1) (Protein p14 TCL1)
Accession No.	Swiss-Prot:P56279NCBI Gene ID:8115
Calculated MW	13 kDa
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1B£B'1500-1B£B'12000

Images



Western blot analysis of K562 mouse-lung lysate, antibody was diluted at 1000. Secondary antibody was diluted at 1:20000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.