

Cytokeratin 4 Polyclonal Antibody

Catalog No: #46782

Package Size: #46782-1 50ul #46782-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Cytokeratin 4 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	The antibody detects endogenous Cytokeratin 4
Immunogen Description	Synthetic peptide from human protein at AA range: 200-260
Conjugates	Unconjugated
Target Name	Cytokeratin?4
Other Names	Keratin, type II cytoskeletal 4 (Cytokeratin-4) (CK-4) (Keratin-4) (K4) (Type-II keratin Kb4)
Accession No.	Swiss-Prot:P19013NCBI Gene ID:3851
Calculated MW	57 kDa
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1:500-1:2000

Images



Western blot analysis of SKOV3 293T lysate, antibody was diluted at 500. Secondary antibody was diluted at 1:20000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.