

RPA14 Polyclonal Antibody

Catalog No: #46786

Package Size: #46786-1 50ul #46786-2 100ul


Description

Product Name	RPA14 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Applications	WB
Species Reactivity	Hu Ms
Specificity	The antibody detects endogenous RPA14
Immunogen Type	protein
Immunogen Description	Synthetic peptide from human protein at AA range: 30-80
Conjugates	Unconjugated
Target Name	RPA14
Other Names	Replication protein A 14 kDa subunit (RP-A p14) (Replication factor A protein 3) (RF-A protein 3)
Accession No.	Swiss-Prot:P35244NCBI Gene ID:6119
Calculated MW	14 kDa
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1B£B1500-1B£B12000

Images



Western blot analysis of lysate, antibody was diluted at 2000.
Secondary antibody was diluted at 1:20000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.