

FSH β Polyclonal Antibody

Catalog No: #46792

Package Size: #46792-1 50ul #46792-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	FSH β Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Applications	IHC
Species Reactivity	Hu Ms Rt
Specificity	The antibody detects endogenous FSH β
Immunogen Type	protein
Immunogen Description	Synthetic peptide from human protein at AA range: 10-50
Conjugates	Unconjugated
Target Name	FSH β
Other Names	Follitropin subunit beta (Follicle-stimulating hormone beta subunit) (FSH-B) (FSH-beta) (Follitropin beta chain)
Accession No.	Swiss-Prot:P01225NCBI Gene ID:2488
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C


Application Details

Immunohistochemistry: IHC 1:50-300

Images



Immunohistochemical analysis of paraffin-embedded human-brain, antibody was diluted at 1:200



Immunohistochemical analysis of paraffin-embedded human-brain, antibody was diluted at 1:200

Note: This product is for in vitro research use only and is not intended for use in humans or animals.