

CD1B Polyclonal Antibody

Catalog No: #46805

Package Size: #46805-1 50ul #46805-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	CD1B Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous CD1B
Immunogen Type	protein
Immunogen Description	Synthetic peptide from human protein at AA range: 60-100
Conjugates	Unconjugated
Target Name	CD1B
Other Names	T-cell surface glycoprotein CD1b (CD antigen CD1b)
Accession No.	Swiss-Prot:P29016NCBI Gene ID:910
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C


Application Details

Immunohistochemistry: IHC 1:50-300


Images



Immunohistochemical analysis of paraffin-embedded human-brain, antibody was diluted at 1:200



Immunohistochemical analysis of paraffin-embedded human-kidney, antibody was diluted at 1:200



Immunohistochemical analysis of paraffin-embedded human-kidney, antibody was diluted at 1:200

Note: This product is for in vitro research use only and is not intended for use in humans or animals.