

CD108 Polyclonal Antibody

Catalog No: #46809

Package Size: #46809-1 50ul #46809-2 100ul


Description

Product Name	CD108 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Applications	IHC
Species Reactivity	Hu Ms Rt
Specificity	The antibody detects endogenous CD108
Immunogen Type	protein
Immunogen Description	Synthetic peptide from human protein at AA range: 101-150
Conjugates	Unconjugated
Target Name	CD108
Other Names	Semaphorin-7A (CDw108) (JMH blood group antigen) (John-Milton-Hargen human blood group Ag) (Semaphorin-K1) (Sema K1) (Semaphorin-L) (Sema L) (CD antigen CD108)
Accession No.	Swiss-Prot:O75326NCBI Gene ID:8482
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Immunohistochemistry: IHC 1:50-300

Images



Immunohistochemical analysis of paraffin-embedded human-brain, antibody was diluted at 1:200

Note: This product is for in vitro research use only and is not intended for use in humans or animals.