

## CD272 Polyclonal Antibody

Catalog No: #46813

Package Size: #46813-1 50ul #46813-2 100ul


## Description

|                       |                                                                                                           |
|-----------------------|-----------------------------------------------------------------------------------------------------------|
| Product Name          | CD272 Polyclonal Antibody                                                                                 |
| Host Species          | Rabbit                                                                                                    |
| Clonality             | Polyclonal                                                                                                |
| Purification          | The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen. |
| Applications          | IHC                                                                                                       |
| Species Reactivity    | Hu                                                                                                        |
| Specificity           | The antibody detects endogenous CD272                                                                     |
| Immunogen Type        | protein                                                                                                   |
| Immunogen Description | Synthetic peptide from human protein at AA range: 101-150                                                 |
| Conjugates            | Unconjugated                                                                                              |
| Target Name           | CD272                                                                                                     |
| Other Names           | B- and T-lymphocyte attenuator (B- and T-lymphocyte-associated protein) (CD antigen CD272)                |
| Accession No.         | Swiss-Prot:Q7Z6A9NCBI Gene ID:151888                                                                      |
| Concentration         | 1 mg/ml                                                                                                   |
| Formulation           | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                   |
| Storage               | Store at -20°C                                                                                            |


## Application Details

Immunohistochemistry: IHC 1:50-300

## Images



Immunohistochemical analysis of paraffin-embedded human-small-intestine, antibody was diluted at 1:200



Immunohistochemical analysis of paraffin-embedded human-small-intestine, antibody was diluted at 1:200

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.