

CKR-1 Polyclonal Antibody

Catalog No: #46822

Package Size: #46822-1 50ul #46822-2 100ul


Description

Product Name	CKR-1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous CKR-1
Immunogen Type	protein
Immunogen Description	Synthetic peptide from human protein at AA range: 151-200
Conjugates	Unconjugated
Target Name	CKR-1
Other Names	C-C chemokine receptor type 1 (C-C CKR-1) (CC-CKR-1) (CCR-1) (CCR1) (HM145) (LD78 receptor) (Macrophage inflammatory protein 1-alpha receptor) (MIP-1alpha-R) (RANTES-R) (CD antigen CD191)
Accession No.	Swiss-Prot:P32246NCBI Gene ID:1230
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C


Application Details

Immunohistochemistry: IHC 1:50-300

Images



Immunohistochemical analysis of paraffin-embedded human-small-intestine, antibody was diluted at 1:200



Immunohistochemical analysis of paraffin-embedded human-spleen, antibody was diluted at 1:200

Note: This product is for in vitro research use only and is not intended for use in humans or animals.