

CD236 Polyclonal Antibody

Catalog No: #46831

Package Size: #46831-1 50ul #46831-2 100ul


Description

Product Name	CD236 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous CD236
Immunogen Type	protein
Immunogen Description	Synthetic peptide from human protein at AA range: 11-60
Conjugates	Unconjugated
Target Name	CD236
Other Names	Glycophorin-C (Glycoconnectin) (Glycophorin-D) (GPD) (Glycoprotein beta) (PAS-2') (Sialoglycoprotein D) (CD antigen CD236)
Accession No.	Swiss-Prot:P04921NCBI Gene ID:2995
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C


Application Details

Immunohistochemistry: IHC 1:50-300


Images



Immunohistochemical analysis of paraffin-embedded human-lung, antibody was diluted at 1:200



Immunohistochemical analysis of paraffin-embedded human-lung, antibody was diluted at 1:200



Immunohistochemical analysis of paraffin-embedded human-spleen, antibody was diluted at 1:200

Note: This product is for in vitro research use only and is not intended for use in humans or animals.