

## CCL14 Polyclonal Antibody

Catalog No: #46832

Package Size: #46832-1 50ul #46832-2 100ul


## Description

|                       |                                                                                                                                                                       |
|-----------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name          | CCL14 Polyclonal Antibody                                                                                                                                             |
| Host Species          | Rabbit                                                                                                                                                                |
| Clonality             | Polyclonal                                                                                                                                                            |
| Purification          | The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.                                                             |
| Applications          | IHC                                                                                                                                                                   |
| Species Reactivity    | Hu                                                                                                                                                                    |
| Specificity           | The antibody detects endogenous CCL14                                                                                                                                 |
| Immunogen Type        | protein                                                                                                                                                               |
| Immunogen Description | Synthetic peptide from human protein at AA range: 44-93                                                                                                               |
| Conjugates            | Unconjugated                                                                                                                                                          |
| Target Name           | CCL14                                                                                                                                                                 |
| Other Names           | C-C motif chemokine 14 (Chemokine CC-1/CC-3) (HCC-1/HCC-3) (HCC-1(1-74)) (NCC-2) (Small-inducible cytokine A14) [Cleaved into: HCC-1(3-74); HCC-1(4-74); HCC-1(9-74)] |
| Accession No.         | Swiss-Prot:Q16627NCBI Gene ID:6358                                                                                                                                    |
| Concentration         | 1 mg/ml                                                                                                                                                               |
| Formulation           | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                                                                               |
| Storage               | Store at -20°C                                                                                                                                                        |


## Application Details

Immunohistochemistry: IHC 1:50-300

## Images



Immunohistochemical analysis of paraffin-embedded human-brain, antibody was diluted at 1:200



Immunohistochemical analysis of paraffin-embedded human-brain, antibody was diluted at 1:200

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.