

IL-27A Polyclonal Antibody

Catalog No: #46853

Package Size: #46853-1 50ul #46853-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	IL-27A Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Applications	IHC
Species Reactivity	Hu Ms Rt
Specificity	The antibody detects endogenous IL-27A
Immunogen Type	protein
Immunogen Description	Synthetic peptide from human protein at AA range: 1-50
Conjugates	Unconjugated
Target Name	IL-27A
Other Names	Interleukin-27 subunit alpha (IL-27 subunit alpha) (IL-27-A) (IL27-A) (p28)
Accession No.	Swiss-Prot:Q8NEV9NCBI Gene ID:246778
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Immunohistochemistry: IHC 1:50-300

Images



Immunohistochemical analysis of paraffin-embedded human-lung, antibody was diluted at 1:200

Note: This product is for in vitro research use only and is not intended for use in humans or animals.