

IL-17B Polyclonal Antibody

Catalog No: #46854

Package Size: #46854-1 50ul #46854-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	IL-17B Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Applications	IHC
Species Reactivity	Hu Rt
Specificity	The antibody detects endogenous IL-17B
Immunogen Type	protein
Immunogen Description	Synthetic peptide from human protein at AA range: 1-50
Conjugates	Unconjugated
Target Name	IL-17B
Other Names	Interleukin-17B (IL-17B) (Cytokine Zcyto7) (Interleukin-20) (IL-20) (Neuronal interleukin-17-related factor)
Accession No.	Swiss-Prot:Q9UHF5NCBI Gene ID:27190
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Immunohistochemistry: IHC 1:50-300

Images



Immunohistochemical analysis of paraffin-embedded human-small-intestine, antibody was diluted at 1:200

Note: This product is for in vitro research use only and is not intended for use in humans or animals.