

IL-17C Polyclonal Antibody

Catalog No: #46856

Package Size: #46856-1 50ul #46856-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	IL-17C Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous IL-17C
Immunogen Type	protein
Immunogen Description	Synthetic peptide from human protein at AA range: 81-130
Conjugates	Unconjugated
Target Name	IL-17C
Other Names	Interleukin-17C (IL-17C) (Cytokine CX2)
Accession No.	Swiss-Prot:Q9P0M4NCBI Gene ID:27189
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C


Application Details

Immunohistochemistry: IHC 1:50-300


Images



Immunohistochemical analysis of paraffin-embedded human-colon, antibody was diluted at 1:200



Immunohistochemical analysis of paraffin-embedded human-colon, antibody was diluted at 1:200



Immunohistochemical analysis of paraffin-embedded human-pancreas, antibody was diluted at 1:200

Note: This product is for in vitro research use only and is not intended for use in humans or animals.