

ZAP-70 Polyclonal Antibody

Catalog No: #46858

Package Size: #46858-1 50ul #46858-2 100ul


Description

Product Name	ZAP-70 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous ZAP-70
Immunogen Type	protein
Immunogen Description	Synthetic peptide from human protein at AA range: 391-440
Conjugates	Unconjugated
Target Name	ZAP-70
Other Names	Tyrosine-protein kinase ZAP-70 (EC 2.7.10.2) (70 kDa zeta-chain associated protein) (Syk-related tyrosine kinase)
Accession No.	Swiss-Prot:P43403NCBI Gene ID:7535
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Immunohistochemistry: IHC 1:50-300

Images



Immunohistochemical analysis of paraffin-embedded human-lymph, antibody was diluted at 1:200

Note: This product is for in vitro research use only and is not intended for use in humans or animals.