

Abin-2 Polyclonal Antibody

Catalog No: #46865

Package Size: #46865-1 50ul #46865-2 100ul


Description

Product Name	Abin-2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Applications	IHC
Species Reactivity	Hu Ms Rt
Specificity	The antibody detects endogenous Abin-2
Immunogen Type	protein
Immunogen Description	Synthetic peptide from human protein at AA range: 61-110
Conjugates	Unconjugated
Target Name	Abin-2
Other Names	TNFAIP3-interacting protein 2 (A20-binding inhibitor of NF-kappa-B activation 2) (ABIN-2) (Fetal liver LKB1-interacting protein)
Accession No.	Swiss-Prot:Q8NFZ5NCBI Gene ID:79155
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C


Application Details

Immunohistochemistry: IHC 1:50-300


Images




Immunohistochemical analysis of paraffin-embedded human-cervical-cancer, antibody was diluted at 1:200



Immunohistochemical analysis of paraffin-embedded human-cervical-cancer, antibody was diluted at 1:200



Immunohistochemical analysis of paraffin-embedded human-liver, antibody was diluted at 1:200



Immunohistochemical analysis of paraffin-embedded human-liver, antibody was diluted at 1:200

Note: This product is for in vitro research use only and is not intended for use in humans or animals.