

TTN Polyclonal Antibody

Catalog No: #46877

Package Size: #46877-1 50ul #46877-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	TTN Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Applications	IHC
Species Reactivity	Hu Ms Rt
Specificity	The antibody detects endogenous TTN
Immunogen Type	protein
Immunogen Description	Synthetic peptide from human protein at AA range: 161-210
Conjugates	Unconjugated
Target Name	TTN
Other Names	Titin (EC 2.7.11.1) (Connectin) (Rhabdomyosarcoma antigen MU-RMS-40.14)
Accession No.	Swiss-Prot:Q8WZ42NCBI Gene ID:7273
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C


Application Details

Immunohistochemistry: IHC 1:50-300


Images




Immunohistochemical analysis of paraffin-embedded human-heart, antibody was diluted at 1:200



Immunohistochemical analysis of paraffin-embedded human-heart, antibody was diluted at 1:200



Immunohistochemical analysis of paraffin-embedded human-skeletal-muscle, antibody was diluted at 1:200



Immunohistochemical analysis of paraffin-embedded human-skeletal-muscle, antibody was diluted at 1:200

Note: This product is for in vitro research use only and is not intended for use in humans or animals.