

ADI1 Polyclonal Antibody

Catalog No: #46879

Package Size: #46879-1 50ul #46879-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ADI1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Applications	IHC
Species Reactivity	Hu Ms Rt
Specificity	The antibody detects endogenous ADI1
Immunogen Type	protein
Immunogen Description	Synthetic peptide from human protein at AA range: 71-120
Conjugates	Unconjugated
Target Name	ADI1
Other Names	1,2-dihydroxy-3-keto-5-methylthiopentene dioxygenase (EC 1.13.11.54) (Acireductone dioxygenase (Fe(2+)-requiring)) (ARD) (Fe-ARD) (Membrane-type 1 matrix metalloproteinase cytoplasmic tail-binding protein 1) (MTCBP-1) (Submergence-induced protein-like factor) (Sip-L)
Accession No.	Swiss-Prot:Q9BV57NCBI Gene ID:55256
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Immunohistochemistry: IHC 1:50-300

Images



Immunohistochemical analysis of paraffin-embedded human-lung, antibody was diluted at 1:200

Note: This product is for in vitro research use only and is not intended for use in humans or animals.