

PF-4 Polyclonal Antibody

Catalog No: #46883

Package Size: #46883-1 50ul #46883-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	PF-4 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Applications	IHC
Species Reactivity	Hu Rt
Specificity	The antibody detects endogenous PF-4
Immunogen Type	protein
Immunogen Description	Synthetic peptide from human protein at AA range: 51-100
Conjugates	Unconjugated
Target Name	PF-4
Other Names	Platelet factor 4 (PF-4) (C-X-C motif chemokine 4) (Iroplact) (Oncostatin-A) [Cleaved into: Platelet factor 4, short form]
Accession No.	Swiss-Prot:P02776NCBI Gene ID:5196
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Immunohistochemistry: IHC 1:50-300

Images



Immunohistochemical analysis of paraffin-embedded human-spleen, antibody was diluted at 1:200

Note: This product is for in vitro research use only and is not intended for use in humans or animals.