

CD22 Polyclonal Antibody

Catalog No: #46886

Package Size: #46886-1 50ul #46886-2 100ul


Description

Product Name	CD22 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu
Specificity	This antibody detects endogenous levels of CD22
Immunogen Type	protein
Immunogen Description	Synthesized peptide derived from human CD22. at AA range: 780-860
Conjugates	Unconjugated
Target Name	CD22
Other Names	B-cell receptor CD22 (B-lymphocyte cell adhesion molecule) (BL-CAM) (Sialic acid-binding Ig-like lectin 2) (Siglec-2) (T-cell surface antigen Leu-14) (CD antigen CD22)
Accession No.	Swiss-Prot:P20273NCBI Gene ID:933
Calculated MW	95 kDa
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1B£B1500-1B£B12000

Images



Western Blot analysis of hela cells using primary antibody diluted at 1:1000(°C overnight). Secondary antibody B£B1Goat Anti-rabbit IgG IRDye 800(diluted at 1:5000, 25B'Γ£C, 1 hour)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.