

CTXL Polyclonal Antibody

Catalog No: #46887

Package Size: #46887-1 50ul #46887-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	CTXL Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu
Specificity	This antibody detects endogenous levels of CTXL
Immunogen Type	protein
Immunogen Description	Synthesized peptide derived from human CTXL. at AA range: 260-290
Conjugates	Unconjugated
Target Name	CTXL
Other Names	V-set and immunoglobulin domain-containing protein 2 (Cortical thymocyte-like protein) (CT-like protein)
Accession No.	Swiss-Prot:Q96IQ7NCBI Gene ID:23584
Calculated MW	34 kDa
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1B£B1500-1B£B12000

Images



Western Blot analysis of hela cells using primary antibody diluted at 1:1000 (°C overnight). Secondary antibody B£B1Goat Anti-rabbit IgG IRDye 800 (diluted at 1:5000, 25B'Γ£C, 1 hour)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.