

AVP Receptor V2 Polyclonal Antibody

Catalog No: #46898

Package Size: #46898-1 50ul #46898-2 100ul



SAB[®]
Signalway Antibody

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	AVP Receptor V2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu Ms
Specificity	AVP Receptor V2 Polyclonal Antibody detects endogenous levels of AVP Receptor V2
Immunogen Type	protein
Immunogen Description	Synthesized peptide derived from human AVP Receptor V2. at AA range: 1-50
Conjugates	Unconjugated
Target Name	AVP Receptor V2
Other Names	Vasopressin V2 receptor (V2R) (AVPR V2) (Antidiuretic hormone receptor) (Renal-type arginine vasopressin receptor)
Accession No.	Swiss-Prot:P30518NCBI Gene ID:554
Calculated MW	57 kDa
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1B£B1500-1B£B12000

Images



Western Blot analysis of 1,mouse-lung 2,mouse-spleen 3,mouse-kidney 4,mouse-heart 5,293 cells using primary antibody diluted at 1:500(°C overnight). Secondary antibodyB£B1Goat Anti-rabbit IgG IRDye 800(diluted at 1:5000, 25B'Γ£C, 1 hour)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.