

FA2H Polyclonal Antibody

Catalog No: #46899

Package Size: #46899-1 50ul #46899-2 100ul



SAB[®]
Signalway Antibody

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	FA2H Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	FA2H Polyclonal Antibody detects endogenous levels of FA2H
Immunogen Type	protein
Immunogen Description	Synthesized peptide derived from human FA2H. at AA range: 101-150
Conjugates	Unconjugated
Target Name	FA2H
Other Names	Fatty acid 2-hydroxylase (EC 1.-.-.-) (Fatty acid alpha-hydroxylase)
Accession No.	Swiss-Prot:Q7L5A8NCBI Gene ID:79152
Calculated MW	55 kDa
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1B£B1500-1B£B12000

Images



Western Blot analysis of mouse-heart cells using primary antibody diluted at 1:2000(°C overnight). Secondary antibodyB£B1Goat Anti-rabbit IgG IRDye 800(diluted at 1:5000, 25B'Γ£C, 1 hour)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.