

Dlx-2 Polyclonal Antibody

Catalog No: #46902

Package Size: #46902-1 50ul #46902-2 100ul



SAB[®]

Signalway Antibody

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	Dlx-2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	Dlx-2 Polyclonal Antibody detects endogenous levels of Dlx-2
Immunogen Type	protein
Immunogen Description	Synthesized peptide derived from human Dlx-2. at AA range: 101-150
Conjugates	Unconjugated
Target Name	Dlx-2
Other Names	Homeobox protein DLX-2
Accession No.	Swiss-Prot:Q07687NCBI Gene ID:1746
Calculated MW	42 kDa
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1B£B1500-1B£B12000

Images



Western Blot analysis of 1,hela 2,mouse-brain cells using primary antibody diluted at 1:500(°C overnight). Secondary antibodyB£B1Goat Anti-rabbit IgG IRDye 800(diluted at 1:5000, 25B'Γ£C, 1 hour)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.