

Jun B Polyclonal Antibody

Catalog No: #46903

Package Size: #46903-1 50ul #46903-2 100ul


Description

|                       |                                                                                                                       |
|-----------------------|-----------------------------------------------------------------------------------------------------------------------|
| Product Name          | Jun B Polyclonal Antibody                                                                                             |
| Host Species          | Rabbit                                                                                                                |
| Clonality             | Polyclonal                                                                                                            |
| Purification          | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Applications          | WB                                                                                                                    |
| Species Reactivity    | Hu Ms Rt                                                                                                              |
| Specificity           | Jun B Polyclonal Antibody detects endogenous levels of Jun B                                                          |
| Immunogen Type        | protein                                                                                                               |
| Immunogen Description | Synthesized peptide derived from human Jun B. at AA range: 231-280                                                    |
| Conjugates            | Unconjugated                                                                                                          |
| Target Name           | Jun B                                                                                                                 |
| Other Names           | Transcription factor jun-B                                                                                            |
| Accession No.         | Swiss-Prot:P17275NCBI Gene ID:3726                                                                                    |
| Calculated MW         | 40 kDa                                                                                                                |
| Concentration         | 1 mg/ml                                                                                                               |
| Formulation           | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                               |
| Storage               | Store at -20°C                                                                                                        |

Application Details

Western blotting: 1B£B1500-1B£B12000

Images



Western Blot analysis of 1,mouse-kidney 2,mouse-heart cells using primary antibody diluted at 1:500(°C overnight). Secondary antibodyB£B1Goat Anti-rabbit IgG IRDye 800( diluted at 1:5000, 25B1°C, 1 hour)

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.