

MYL2 Polyclonal Antibody

Catalog No: #46911

Package Size: #46911-1 50ul #46911-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	MYL2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	MYL2 Polyclonal Antibody detects endogenous levels of MYL2
Immunogen Type	protein
Immunogen Description	Synthesized peptide derived from human MYL2. at AA range: 91-140
Conjugates	Unconjugated
Target Name	MYL2
Other Names	Myosin regulatory light chain 2, ventricular/cardiac muscle isoform (MLC-2) (MLC-2v)
Accession No.	Swiss-Prot:P10916NCBI Gene ID:4633
Calculated MW	18 kDa
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1B£B1500-1B£B12000

Images



Western Blot analysis of 1,mouse-heart 2,Hela cells using primary antibody diluted at 1:500(°C overnight). Secondary antibodyB£B1Goat Anti-rabbit IgG IRDye 800(diluted at 1:5000, 25B'Γ£C, 1 hour)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.