

Cytochrome b Polyclonal Antibody

Catalog No: #46912

Package Size: #46912-1 50ul #46912-2 100ul


Description

Product Name	Cytochrome b Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu Ms
Specificity	Cytochrome b Polyclonal Antibody detects endogenous levels of Cytochrome b
Immunogen Type	protein
Immunogen Description	Synthesized peptide derived from human Cytochrome b. at AA range: 331-380
Conjugates	Unconjugated
Target Name	Cytochrome b
Other Names	Cytochrome b (Complex III subunit 3) (Complex III subunit III) (Cytochrome b-c1 complex subunit 3) (Ubiquinol-cytochrome-c reductase complex cytochrome b subunit)
Accession No.	Swiss-Prot:P00156NCBI Gene ID:4519
Calculated MW	48 kDa
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1B£B'1500-1B£B'12000

Images



Western Blot analysis of 1,mouse-lung 2,mouse-brain 3,mouse-spleen 4,mouse-kidney 5,mouse-heart cells using primary antibody diluted at 1:500(°C overnight). Secondary antibodyB£B'Goat Anti-rabbit IgG IRDye 800(diluted at 1:5000, 25B'Γ£C, 1 hour)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.