

C8 α Polyclonal Antibody

Catalog No: #46913

Package Size: #46913-1 50ul #46913-2 100ul



SAB[®]

Signalway Antibody

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	C8 α Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu
Specificity	C8 α Polyclonal Antibody detects endogenous levels of C8 α
Immunogen Type	protein
Immunogen Description	Synthesized peptide derived from human C8 α. at AA range: 201-250
Conjugates	Unconjugated
Target Name	C8 α
Other Names	Complement component C8 alpha chain (Complement component 8 subunit alpha)
Accession No.	Swiss-Prot:P07357NCBI Gene ID:731
Calculated MW	60 kDa
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1B£B1500-1B£B12000

Images



Western Blot analysis of 1,mouse-kidney 2,mouse-heart cells using primary antibody diluted at 1:500(°C overnight). Secondary antibodyB£B1Goat Anti-rabbit IgG IRDye 800(diluted at 1:5000, 25B1°C, 1 hour)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.