

## AVPR1A Polyclonal Antibody

Catalog No: #46914

Package Size: #46914-1 50ul #46914-2 100ul


## Description

|                       |                                                                                                                                     |
|-----------------------|-------------------------------------------------------------------------------------------------------------------------------------|
| Product Name          | AVPR1A Polyclonal Antibody                                                                                                          |
| Host Species          | Rabbit                                                                                                                              |
| Clonality             | Polyclonal                                                                                                                          |
| Purification          | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.               |
| Applications          | WB                                                                                                                                  |
| Species Reactivity    | Hu Ms                                                                                                                               |
| Specificity           | AVPR1A Polyclonal Antibody detects endogenous levels of AVPR1A                                                                      |
| Immunogen Type        | protein                                                                                                                             |
| Immunogen Description | Synthesized peptide derived from human AVPR1A. at AA range: 11-60                                                                   |
| Conjugates            | Unconjugated                                                                                                                        |
| Target Name           | AVPR1A                                                                                                                              |
| Other Names           | Vasopressin V1a receptor (V1aR) (AVPR V1a) (Antidiuretic hormone receptor 1a) (Vascular/hepatic-type arginine vasopressin receptor) |
| Accession No.         | Swiss-Prot:P37288NCBI Gene ID:552                                                                                                   |
| Calculated MW         | 55 kDa                                                                                                                              |
| Concentration         | 1 mg/ml                                                                                                                             |
| Formulation           | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                                             |
| Storage               | Store at -20°C                                                                                                                      |

## Application Details

Western blotting: 1B£B'1500-1B£B'12000

## Images



Western Blot analysis of mouse-lung cells using primary antibody diluted at 1:500(°C overnight). Secondary antibodyB£B'Goat Anti-rabbit IgG IRDye 800( diluted at 1:5000, 25B'£C, 1 hour)

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.