

# AFG3L2 Polyclonal Antibody

Catalog No: #46915

Package Size: #46915-1 50ul #46915-2 100ul


## Description

|                       |                                                                                                                       |
|-----------------------|-----------------------------------------------------------------------------------------------------------------------|
| Product Name          | AFG3L2 Polyclonal Antibody                                                                                            |
| Host Species          | Rabbit                                                                                                                |
| Clonality             | Polyclonal                                                                                                            |
| Purification          | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Applications          | WB                                                                                                                    |
| Species Reactivity    | Hu                                                                                                                    |
| Specificity           | AFG3L2 Polyclonal Antibody detects endogenous levels of AFG3L2                                                        |
| Immunogen Type        | protein                                                                                                               |
| Immunogen Description | Synthesized peptide derived from human AFG3L2. at AA range: 744-793                                                   |
| Conjugates            | Unconjugated                                                                                                          |
| Target Name           | AFG3L2                                                                                                                |
| Other Names           | AFG3-like protein 2 (EC 3.4.24.-) (Paraplegin-like protein)                                                           |
| Accession No.         | Swiss-Prot:Q9Y4W6NCBI Gene ID:10939                                                                                   |
| Calculated MW         | 88 kDa                                                                                                                |
| Concentration         | 1 mg/ml                                                                                                               |
| Formulation           | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                               |
| Storage               | Store at -20°C                                                                                                        |

## Application Details

Western blotting: 1B£B'1500-1B£B'12000

## Images



Western Blot analysis of 1,mouse-heart 2,293T cells using primary antibody diluted at 1:500(°C overnight). Secondary antibodyB£B'1Goat Anti-rabbit IgG IRDye 800( diluted at 1:5000, 25B'Γ£C, 1 hour)

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.