

NMUR2 Polyclonal Antibody

Catalog No: #46916

Package Size: #46916-1 50ul #46916-2 100ul



SAB[®]

Signalway Antibody

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	NMUR2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu
Specificity	NMUR2 Polyclonal Antibody detects endogenous levels of NMUR2
Immunogen Type	protein
Immunogen Description	Synthesized peptide derived from human NMUR2. at AA range: 1-50
Conjugates	Unconjugated
Target Name	NMUR2
Other Names	Neuromedin-U receptor 2 (NMU-R2) (G-protein coupled receptor FM-4) (G-protein coupled receptor TGR-1)
Accession No.	Swiss-Prot:Q9GZQ4NCBI Gene ID:56923
Calculated MW	46 kDa
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1B£B1500-1B£B12000

Images



Western Blot analysis of mouse-brain cells using primary antibody diluted at 1:2000(°C overnight). Secondary antibodyB£B1Goat Anti-rabbit IgG IRDye 800(diluted at 1:5000, 25B'Γ£C, 1 hour)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.