

Maspin Polyclonal Antibody

Catalog No: #46917

Package Size: #46917-1 50ul #46917-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Maspin Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	Maspin Polyclonal Antibody detects endogenous levels of Maspin
Immunogen Type	protein
Immunogen Description	Synthesized peptide derived from human Maspin. at AA range: 201-250
Conjugates	Unconjugated
Target Name	Maspin
Other Names	Serpin B5 (Maspin) (Peptidase inhibitor 5) (PI-5)
Accession No.	Swiss-Prot:P36952NCBI Gene ID:5268
Calculated MW	42 kDa
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1B£B1500-1B£B12000

Images



Western Blot analysis of 293T cells using primary antibody diluted at 1:2000(°C overnight). Secondary antibodyB£B1Goat Anti-rabbit IgG IRDye 800(diluted at 1:5000, 25B'£C, 1 hour)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.