

COL13A1 Polyclonal Antibody

Catalog No: #46920

Package Size: #46920-1 50ul #46920-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	COL13A1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu Ms
Specificity	COL13A1 Polyclonal Antibody detects endogenous levels of COL13A1
Immunogen Type	protein
Immunogen Description	Synthesized peptide derived from human COL13A1. at AA range: 41-90
Conjugates	Unconjugated
Target Name	COL13A1
Other Names	Collagen alpha-1(XIII) chain (COLXIII A1)
Accession No.	Swiss-Prot:Q5TAT6NCBI Gene ID:1305
Calculated MW	70 kDa
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1B£B1500-1B£B12000

Images



Western Blot analysis of 1,Hela 2,Mouse-kidney cells using primary antibody diluted at 1:2000(°C overnight). Secondary antibodyB£B1Goat Anti-rabbit IgG IRDye 800(diluted at 1:5000, 25B'Γ£C, 1 hour)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.