

C1RL Polyclonal Antibody

Catalog No: #46921

Package Size: #46921-1 50ul #46921-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	C1RL Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu
Specificity	C1RL Polyclonal Antibody detects endogenous levels of C1RL
Immunogen Type	protein
Immunogen Description	Synthesized peptide derived from human C1RL. at AA range: 371-420
Conjugates	Unconjugated
Target Name	C1RL
Other Names	Complement C1r subcomponent-like protein (C1r-LP) (C1r-like protein) (EC 3.4.21.-) (C1r-like serine protease analog protein) (CLSPa)
Accession No.	Swiss-Prot:Q9NZP8NCBI Gene ID:51279
Calculated MW	53 kDa
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1B£B'1500-1B£B'12000

Images



Western Blot analysis of hela cells using primary antibody diluted at 1:500(°C overnight). Secondary antibodyB£B'Goat Anti-rabbit IgG IRDye 800(diluted at 1:5000, 25B'Γ£C, 1 hour)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.