

## PTGR2 Polyclonal Antibody

Catalog No: #47347

Package Size: #47347-1 50ul #47347-2 100ul


## Description

|                       |                                                                                                                       |
|-----------------------|-----------------------------------------------------------------------------------------------------------------------|
| Product Name          | PTGR2 Polyclonal Antibody                                                                                             |
| Host Species          | Rabbit                                                                                                                |
| Clonality             | Polyclonal                                                                                                            |
| Isotype               | IgG                                                                                                                   |
| Purification          | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Applications          | WB;IHC;IF;ELISA                                                                                                       |
| Species Reactivity    | Human;Rat;Mouse                                                                                                       |
| Specificity           | PTGR2 Polyclonal Antibody detects endogenous levels of PTGR2 protein.                                                 |
| Immunogen Description | The antiserum was produced against synthesized peptide derived from human ZADH1.                                      |
| Conjugates            | Unconjugated                                                                                                          |
| Target Name           | PTGR2                                                                                                                 |
| Calculated MW         | 35kD                                                                                                                  |
| Concentration         | 1 mg/ml                                                                                                               |
| Formulation           | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                               |
| Storage               | Store at -20°C/1 year                                                                                                 |


## Application Details

WB: 1/500 - 1/2000. IHC: 1/100 - 1/300. ELISA: 1/10000.. IF 1:50-200

## Images



Western Blot analysis of K562 cells using PTGR2 Polyclonal Antibody



Immunohistochemical analysis of paraffin-embedded human tonsil.

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.